


From Eye to Insight



The Perfect Aid for Placement of Premium IOLs

Leica ToricEyePiece



POSITION TORIC IOLS QUICKLY AND EASILY

Cost-effective, easy-to-use aid for time savings in cataract surgery of premium IOLs



PERFECT AID FOR IOL POSITIONING

The rotatable graticule* is superimposed over the microscope image and helps the surgeon find the correct IOL positioning angle.

EASY TO USE

Simply rotate the graticule using the sterile handle and identify the desired target position for the toric IOL.

TIME-SAVING

Typically two ink markings are necessary to implant a toric IOL: one ink marking before surgery to mark the zero axis and another during surgery to find the target axis for implantation. With the Leica ToricEyePiece, the second marking may not be necessary, saving the surgeon time.

COST-EFFECTIVE

In comparison to expensive software systems, the Leica ToricEyePiece is the ideal, cost-effective solution for cataract surgeons.



TECHNICAL DATA

Eyepiece	10× / 21B, Type II
Sterility	Autoclavable components: sterile lever
Compatibility	All Leica ophthalmic microscopes with 10× / 21B, Type II eyepieces

*graticule not calibrated



Leica Microsystems (Schweiz) AG
Max Schmidheiny-Strasse 201
9435 Heerbrugg, Switzerland



Leica ToricEyePiece is a Class I medical device.

Not all products or services are approved or offered in every market. Approved labeling and instructions may vary between countries. Please contact your local Leica representative for details.

CONNECT
WITH US!



Leica Microsystems (Schweiz) AG · Max Schmidheiny Strasse 201 · CH-9435 Heerbrugg

T +41 71 726 3333 · F +41 71 726 3399

www.leica-microsystems.com